

PROCESS AND STERILE FILTRATION



AIR AND GAS TREATMENT

Since 1988 OMEGA AIR is developing and manufacturing quality solutions for treatment of industrial as well as process air and gas. Due to constant investment in research and development and in new production technologies we can offer one of the most comprehensive ranges of products for general purpose process filtration, steam filtration, sterile filtration and sterile ventilation of tanks.

TYPICAL APPLICATIONS


- Aseptic processing
- Aseptic packaging
- PET bottle blowing
- Moving ingredients
- Tank ventilation
- General filtration applications
-

TYPICAL INDUSTRIES

- Food and Beverage
- Pharmaceutical / Biopharmaceutical
- Electronics / Semiconductor
- Medical / Health care
-

CONVENTIONAL COMPRESSED AIR SYSTEM

In process industry purity requirements for air/gas are usually higher or at least different compared to purity requirements for standard industrial applications. Production of compressed air in process industry usually starts in conventional compressed air system with standard equipment for treatment of industrial compressed air/gas. At certain point there is a switchover from industrial to process filtration and from this point on different materials and filtration procedures are applied.




PROCESS AIR/GAS SYSTEM

Due to contaminants such as solid particles, moisture, hydrocarbons or bacteria compressed air/gas when used in process can harm and cause quality and efficiency problems. Applications in process industry require compressed air or gas where above mentioned contaminants are removed or reduced to acceptable level for specific application.

LEGEND

- 1 Stainless steel process filter, PN filter element
- 2 Stainless steel process filter, PPF filter element
- 3 Stainless steel sterile filter, VSF filter element
- 4 Steam generator
- 5 Steam filter, PIW filter element




VENTILATION OF TANKS

During filling and emptying of tank with liquid ambient air is being vented in and out of the tank. In order not to harm the process by letting contaminants from ambient into the tank “air-vent” filters (AV series) are installed on the top of the tank. Inside is filter element (AVF) which is designed to intercept any potential particulate or aerosol contaminants including bacteria.

LEGEND

- 1 Tank (pressure vessel)
- 2 AV filter




WHAT IS STERILE AIR

In process filtration applications “sterile” refers to “free from live bacteria or other microorganisms”.


SIZE OF BACTERIA

Size Comparison Chart




HOW TO PRODUCE STERILE AIR

To produce sterile air proper sterile filter needs to be installed (e.g. SF filter housing with VSF filter element) to filter out / intercept bacteria. It is essential that air on the inlet side of sterile filter is adequately prefiltered upstream and that complete “sterile filter installation” is sterilized regularly. Sterilization is required because intercepted bacteria can colonize in the surface of filter media and begin to grow/migrate through filter media. Two most common types of sterilization are:



STERILIZATION IN PLACE (SIP) where filter or even larger segments of installation are sterilized with stream of steam flowing through complete filter installation. This procedure is done in place of use without removing any parts from the installation.

AUTOClave sterilization where filter element or in some cases even complete filter housing is removed from installation and sterilized in autoclave. Autoclave is a device used to sterilize products by subjecting them to high pressure saturated steam.




Frequency of sterilization depend on specific application requirements.

STERILIZATION IN PLACE PROCEDURE


Usually every process or application has specific operating requirements and also sterilization is then performed in accordance with this specific requirements.

Cumulative steaming time:


- 121 °C/250 °F, Sterilization 30 min, Heating and cooling 30 min (100 cycles)
- 132 °C/270 °F, Sterilization 20 min, Heating and cooling 40 min (100 cycles)
- 143 °C/290 °F, Sterilization 10 min, Heating and cooling 50 min (100 cycles)




Sterilization cycle



Steam Saturation Curve



Forward steam in place



Reverse steam in place



14 bar
operating pressure

75 to 21120 Nm³/h
volume flow rate

1/4" to DN200
connections

up to 150 °C
operating temperature range

stainless steel 1.4404-standard
stainless steel 1.4301-option
material

DESCRIPTION

PF process filters are designed for applications in process industry, where the risk for corrosion of compressed air system components is very high.

Required compressed air quality according to standard ISO 8571-1 can be achieved with 9 different grades of filter elements.

PF process filter housing can be used in variety of applications. For applications not listed please contact producer or your local distributor.

For any other technical gas please contact producer or your local distributor.

For oil removal, coalescing filter element must be installed and flow direction inside-out must be provided. General arrangement is bowl on the top and filter head on the bottom.


Fluid group 1 on request.

APPLICATIONS


- Packing industry
- Biotechnology
- Breweries
- Chemical industry
- Dairies
- Fermentation processes
- Food and beverage industry
- Pharmaceutical industry
- Hospitals
- Steam

PF SERIES

STAINLESS STEEL PROCESS FILTERS





TECHNICAL DATA										FILTER ELEMENTS									
Filter housing size	Pipe size D	Oger. press.	Flow rate at 7 bar(g), 20 °C		Dimensions [mm]				Mass	PI prefilter 1; 20 µm	PIW prefilter 1; 20 µm	PN prefilter 5; 25 µm	PP prefilter 3 µm	PR prefilter 1 µm	PM microfilter 0,1 µm	PS microfilter 0,01 µm	PA activated carbon	CKL-P	
	inch		bar/psi	Nm³/h	scfm	A	B	C											E
PF 005	1/4"	14/203	75	44	206	120	76,1	1/4"	1,8	0310 PI	0310 PIW	0310 PN	0310 PP	0310 PR	0310 PM	0310 PS	0310 PA	0310 PA	
PF 007	3/8"	14/203	105	62	236	120	76,1	1/4"	2,0	0410 PI	0410 PIW	0410 PN	0410 PP	0410 PR	0410 PM	0410 PS	0410 PA	0410 PA	
PF 010	1/2"	14/203	150	88	239	121	76,1	1/4"	2,1	0420 PI	0420 PIW	0420 PN	0420 PP	0420 PR	0420 PM	0420 PS	0420 PA	0420 PA	
PF 018	3/4"	14/203	225	132	263	121	76,1	1/4"	2,2	0520 PI	0520 PIW	0520 PN	0520 PP	0520 PR	0520 PM	0520 PS	0520 PA	0520 PA	
PF 030	1"	14/203	315	185	280	136	88,9	1/4"	3,0	0525 PI	0525 PIW	0525 PN	0525 PP	0525 PR	0525 PM	0525 PS	0525 PA	0525 PA	
PF 047	1 1/4"	14/203	420	247	343	155	88,9	1/4"	3,4	0725 PI	0725 PIW	0725 PN	0725 PP	0725 PR	0725 PM	0725 PS	0725 PA	0725 PA	
PF 070	1 1/2"	14/203	600	353	376	180	114,3	1/4"	4,6	0730 PI	0730 PIW	0730 PN	0730 PP	0730 PR	0730 PM	0730 PS	0730 PA	0730 PA	
PF 094	2"	14/203	900	530	445	180	114,3	1/4"	5,2	1030 PI	1030 PIW	1030 PN	1030 PP	1030 PR	1030 PM	1030 PS	1030 PA	1030 PA	
PF 150	2"	14/203	1.260	742	572	180	114,3	1/4"	6,0	1530 PI	1530 PIW	1530 PN	1530 PP	1530 PR	1530 PM	1530 PS	1530 PA	1530 PA	
PF 175	2 1/2"	14/203	1.680	989	736	226	139,7	1/4"	9,6	2030 PI	2030 PIW	2030 PN	2030 PP	2030 PR	2030 PM	2030 PS	2030 PA	2030 PA	
PF 200	3"	14/203	2.400	1.413	979	226	139,7	1/4"	13,7	3030 PI	3030 PIW	3030 PN	3030 PP	3030 PR	3030 PM	3030 PS	3030 PA	3030 PA	
PF 240	3"	14/203	3.600	2.119	1041	256	168,3	1/4"	18,5	3050 PI	3050 PIW	3050 PN	3050 PP	3050 PR	3050 PM	3050 PS	3050 PA	3050 PA	
PF 450	DN100	10/145	5.040	2.966	981	410	219,1	1"	56	3x2030 PI	3x2030 PIW	3x2030 PN	3x2030 PP	3x2030 PR	3x2030 PM	3x2030 PS	3x2030 PA	-	
PF 600	DN100	10/145	6.720	3.955	1240	410	219,1	1"	60	3x3030 PI	3x3030 PIW	3x3030 PN	3x3030 PP	3x3030 PR	3x3030 PM	3x3030 PS	3x3030 PA	-	
PF 900	DN150	10/145	9.600	5.650	1311	480	273,0	1"	90	4x3030 PI	4x3030 PIW	4x3030 PN	4x3030 PP	4x3030 PR	4x3030 PM	4x3030 PS	4x3030 PA	-	
PF 1200	DN150	10/145	13.440	7.910	1330	540	323,9	1"	112	6x3030 PI	6x3030 PIW	6x3030 PN	6x3030 PP	6x3030 PR	6x3030 PM	6x3030 PS	6x3030 PA	-	
PF 1800	DN200	10/145	17.280	10.171	1496	660	406,4	1"	201	8x3030 PI	8x3030 PIW	8x3030 PN	8x3030 PP	8x3030 PR	8x3030 PM	8x3030 PS	8x3030 PA	-	
PF 2000	DN200	10/145	21.120	12.431	1496	660	406,4	1"	202	10x3030 PI	10x3030 PIW	10x3030 PN	10x3030 PP	10x3030 PR	10x3030 PM	10x3030 PS	10x3030 PA	-	
										quality class - solids (ISO 8573-1)	-	-	-	6	3	2	1	1 ¹⁾	-
										quality class - oils (ISO 8573-1)	-	-	-	-	-	2	1	1	-
										pressure drop - new element-dry [mbar]	≤2600; ≤60	≤2600; ≤60	10	10	20	50	80	60	-
										filter media	sintered INOX 1.4404	sintered INOX 1.4404	stainless steel mesh 1.4301	acrylic fibres, cellulose		borosilicate micro fibres		borosilicate micro fibres, activ. carbon	-
										pleated version	-	-	-	✓	✓	✓	✓	-	-
										wrapped version	-	-	✓	-	-	-	-	✓	-
										sintered version	✓	✓	-	-	-	-	-	-	-
min. operating temperature (°C / °F)										1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35
max. operating temperature (°C / °F)										150 / 302	150 / 302	150 / 302	65 / 149	120 / 248	120 / 248	120 / 248	45 / 113	120 / 248	

CORRECTION FACTORS

Operating pressure [bar]	2	3	4	5	6	7	8	9	10	11	12	13	14
Operating pressure [psi]	29	44	58	72	87	100	115	130	145	160	174	189	203
Correction factor	0,38	0,50	0,63	0,75	0,88	1	1,13	1,25	1,38	1,50	1,63	1,75	1,88

¹⁾ Valid if "S" filter cartridge is installed upstream.



50 bar
operating pressure

150 to 2400 Nm³/h
volume flow rate

1/2" to 3"
connections

1,5 to 150 °C
operating temperature range

stainless steel 1.4404-standard
stainless steel 1.4301-option
material

DESCRIPTION

HPF process filters are designed for high pressure applications in process industry, where the risk for corrosion of compressed air system components is very high.

Required compressed air quality according to standard ISO 8571-1 can be achieved with 9 different grades of filter elements.

HPF process filter housing can be used in variety of applications. For applications not listed please contact producer or your local distributor.

For any other technical gas please contact producer or your local distributor.
For oil removal, coalescing filter element must be installed and flow direction inside-out must be provided. General arrangement is bowl on the top and filter head on the bottom.


Fluid group 1 on request.

APPLICATIONS

- Packing industry
- Biotechnology
- Breweries
- Chemical industry
- Dairies
- Fermentation processes
- Food and beverage industry
- Pharmaceutical industry
- Hospitals

HPF SERIES

HIGH PRESSURE STAINLESS STEEL PROCESS FILTERS




AOK 50SS



TD 420M





TECHNICAL DATA									FILTER ELEMENTS									
Filter housing size	Pipe size D	Oper. press.	Flow rate at 7 bar(g), 20 °C		Dimensions [mm]			Mass	PI prefilter 1; 20 µm	PIW prefilter 1; 20 µm	PN prefilter 5; 25 µm	PP prefilter 3 µm	PR prefilter 1 µm	PM microfilter 0,1 µm	P5 microfilter 0,01 µm	PA activated carbon	CKL-P	
	inch								bar/psi	Nm³/h	scfm	A	B	C	kg			
HPF 010	1/2"	50/725	150	88	243	121	76,1	2,5	0420 PI	0420 PIW	0420 PN	0420 PP	0420 PR	0420 PM	0420 PS	0420 PA	0420 CKL-P	
HPF 018	3/4"	50/725	225	132	267	129	76,1	2,7	0520 PI	0520 PIW	0520 PN	0520 PP	0520 PR	0520 PM	0520 PS	0520 PA	0520 CKL-P	
HPF 030	1"	50/725	315	185	286	143	88,9	3,4	0525 PI	0525 PIW	0525 PN	0525 PP	0525 PR	0525 PM	0525 PS	0525 PA	0525 CKL-P	
HPF 047	1 1/4"	50/725	420	247	336	160	88,9	3,9	0725 PI	0725 PIW	0725 PN	0725 PP	0725 PR	0725 PM	0725 PS	0725 PA	0725 CKL-P	
HPF 070	1 1/2"	50/725	600	353	400	182	114,3	5,6	0730 PI	0730 PIW	0730 PN	0730 PP	0730 PR	0730 PM	0730 PS	0730 PA	0730 CKL-P	
HPF 094	2"	50/725	900	530	481	180	114,3	9,2	1030 PI	1030 PIW	1030 PN	1030 PP	1030 PR	1030 PM	1030 PS	1030 PA	1030 CKL-P	
HPF 150	2"	50/725	1260	742	607	180	114,3	10,9	1530 PI	1530 PIW	1530 PN	1530 PP	1530 PR	1530 PM	1530 PS	1530 PA	1530 CKL-P	
HPF 200	3"	50/725	2400	1413	1018	228	139,7	11,5	3030 PI	3030 PIW	3030 PN	3030 PP	3030 PR	3030 PM	3030 PS	3030 PA	3030 CKL-P	
	quality class - solids (ISO 8573-1)								-	-	-	6	3	2	1	1 ¹⁾	-	
	quality class - oils (ISO 8573-1)								-	-	-	-	-	2	1	1	-	
	pressure drop - new element-dry [mbar]								≤2600; ≤60	≤2600; ≤60	10	10	20	50	80	60	-	
	filter media								sintered INOX 1.4404	sintered INOX 1.4404	stainless steel mesh 1.4301	acrylic fibres, cellulose		borosilicate micro fibres		borosilicate micro fibres, activ. carbon	-	
	pleated version								-	-	-	✓	✓	✓	✓	-	-	
	wrapped version								-	-	✓	-	-	-	-	✓	-	
	sintered version								✓	✓	-	-	-	-	-	-	-	
	min. operating temperature (°C / °F)								1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35
	max. operating temperature (°C / °F)								150 / 302	150 / 302	150 / 302	65 / 149	120 / 248	120 / 248	120 / 248	45 / 113	120 / 248	

CORRECTION FACTORS									
Operating pressure [bar]	3	5	7	10	13	16	20	30	40
Operating pressure [psi]	44	72	100	145	189	232	290	435	580
Correction factor	0,50	0,75	1	1,38	1,75	2,13	2,63	3,88	5,13

¹⁾ Valid if "S" filter cartridge is installed upstream.



10 bar
operating pressure

3 to 600 l/min
volume flow rate

3/8" to 3"
connections

1,5 to 150 °C
operating temperature range

stainless steel 1.4404
material

DESCRIPTION

PFL process filter housings have been specifically developed for liquid filtration applications.

To meet the required liquid quality, appropriate filter element must be installed into filter housing.

For any other technical gas/liquid, please contact us or your local dealer.

PFL process filter housing can be used in variety of applications. For applications not listed please contact producer or your local dealer.

APPLICATIONS

- Biotechnology
- Breweries
- Chemical industry
- Dairies
- Fermentation processes
- Food and beverage industry
- Hospitals

PFL SERIES

STAINLESS STEEL PROCESS FILTERS FOR LIQUIDS





TECHNICAL DATA

Filter housing size	Pipe size	Operating pressure	Flow rate at 7 bar(g), 20 °C		Dimensions [mm]				Mass
	inch		l/min	gpm	A	B	C	D	kg
PFL 003	3/8"	10/145	3	0,8	270	110	76,1	90	1,2
PFL 008	3/8"	10/145	8	2,1	330	110	76,1	90	1,4
PFL 012	1"	10/145	12	3,2	402	220	114,3	105	4,4
PFL 025	1"	10/145	25	6,6	540	220	114,3	105	5,1
PFL 050	1"	10/145	50	13,2	785	220	114,3	105	6,4
PFL 075	1"	10/145	75	19,8	1.029	220	114,3	105	7,8
PFL 080	1 1/2"	10/145	80	21,1	680	320	168,3	192	16
PFL 150	1 1/2"	10/145	150	39,6	930	320	168,3	192	18
PFL 225	1 1/2"	10/145	225	59,4	1.200	320	168,3	192	20
PFL 250	2"	10/145	250	66	960	400	219,1	165	24
PFL 375	2"	10/145	375	99	1.210	400	219,1	165	27
PFL 400*	2 1/2"	10/145	400	106	980	500	323,9	175	36
PFL 600*	3"	10/145	600	159	1.230	500	323,9	175	41


CORRECTION FACTORS

Operating pressure [bar]	2	3	4	5	6	7	8	9	10
Operating pressure [psi]	29	44	58	72	87	100	115	130	145
Correction factor	0,38	0,50	0,63	0,75	0,88	1	1,13	1,25	1,38

Flow capacity based on water.

Standard is BSP pipe connection, other pipe connection on request.

* Max. operating pressure for liquids is 7 bar; 3 bar for saturated steam.





10 bar
operating pressure

3 to 600 l/min
volume flow rate

DN10 to DN65
connections

1,5 to 150 °C
operating temperature range

stainless steel 1.4404
material

DESCRIPTION

SFL stainless steel sterile filter housings have been specifically developed for liquid filtration applications.

To meet the required liquid purity, appropriate filter element must be installed into filter housing.

For any other technical gas/liquid, please contact us or your local dealer.

SFL sterile filter housing can be used in variety of applications. For applications not listed please contact us or your local dealer.

APPLICATIONS

- Biotechnology
- Breweries
- Chemical industry
- Dairies
- Fermentation processes
- Food and beverage industry
- Hospitals

SFL SERIES

STAINLESS STEEL STERILE FILTERS FOR LIQUIDS






TECHNICAL DATA									
Filter housing size	Pipe size	Operating pressure	Flow rate at 7 bar(g), 20 °C		Dimensions [mm]				Mass
	inch	bar/psi	l/min	gpm	A	B	C	D	kg
SFL 003	DN10	10/145	3	0,8	260	110	76,1	90	1,2
SFL 008	DN10	10/145	8	2,1	323	110	76,1	90	1,4
SFL 012	DN25	10/145	12	3,2	396	220	114,3	105	4,4
SFL 025	DN25	10/145	25	6,6	530	220	114,3	105	5,1
SFL 050	DN25	10/145	50	13,2	775	220	114,3	105	6,4
SFL 075	DN25	10/145	75	19,8	1.019	220	114,3	105	7,8
SFL 080	DN40	10/145	80	21,1	680	320	168,3	192	16
SFL 150	DN40	10/145	150	39,6	930	320	168,3	192	18
SFL 225	DN40	10/145	225	59,4	1.200	320	168,3	192	20
SFL 250	DN50	10/145	250	66	960	400	219,1	165	24
SFL 375	DN50	10/145	375	99	1.210	400	219,1	165	27
SFL 400*	DN65	10/145	400	106	980	500	323,9	175	36
SFL 600*	DN65	10/145	600	159	1.230	500	323,9	175	41

CORRECTION FACTORS									
Operating pressure [bar]	2	3	4	5	6	7	8	9	10
Operating pressure [psi]	29	44	58	72	87	100	115	130	145
Correction factor	0,38	0,50	0,63	0,75	0,88	1	1,13	1,25	1,38

Flow capacity based on water.

Standard is BSP pipe connection, other pipe connection on request.

* Max. operating pressure for liquids is 7 bar; 3 bar for saturated steam.





10 bar
operating pressure

75 to 21.120 Nm³/h
volume flow rate

DN10 to DN80 TC ISO
DN100 to DN200 EN
connections

1,5°C to 150°C
operating temperature range

stainless steel 1.4404-standard
stainless steel 1.4301-option
material

DESCRIPTION

SF stainless steel sterile filters are designed for high-efficient sterile filtration of compressed air, process air and technical gasses. SF sterile filters provide high-efficient removal of sub-micron particles down to 0,01µm including microorganisms (bacteria). Filter housing surface is polished down to grade Ra0,8.

Required compressed air quality according to standard ISO 8571-1 can be achieved with appropriate filter element.

For any other technical gas please contact producer or your local distributor.

Fluid group 1 on request.

APPLICATIONS


- Packing industry
- Biotechnology
- Breweries
- Chemical industry
- Dairies
- Fermentation processes
- Food and beverage industry
- Pharmaceutical industry
- Hospitals

SF SERIES

STAINLESS STEEL STERILE FILTERS





TECHNICAL DATA											FILTER ELEMENTS	
Filter housing size	Pipe size	Max. oper. pressure	Flow rate at 7 bar(g), 20 °C		Dimensions [mm]					Mass	VSF sterile 0,01 µm	
	inch	bar/psi	Nm³/h	scfm	A	B	C	D	E	kg		
SF 0006	DN10	10/145	75	44	220	130	76,1	74	1/8"	1,6	1 x 0310 VSF	
SF 0009	DN10	10/145	105	62	248	130	76,1	74	1/8"	1,6	1 x 0410 VSF	
SF 0012	DN15	10/145	150	88	256	128	76,1	78	1/8"	1,8	1 x 0420 VSF	
SF 0018	DN15	10/145	225	132	288	128	76,1	78	1/8"	1,8	1 x 0520 VSF	
SF 0032	DN25	10/145	315	185	315	180	114,3	99	1/4"	3,1	1 x 0530 VSF	
SF 0048	DN32	10/145	600	353	368	181	114,3	90	1/4"	3,4	1 x 0730 VSF	
SF 0072	DN40	10/145	900	530	449	174	114,3	65	1/4"	3,6	1 x 1030 VSF	
SF 0108	DN50	10/145	1.260	742	591	188	114,3	110	1/4"	4,9	1 x 1530 VSF	
SF 0144	DN65	10/145	1.680	989	741	205	139,7	136	1/4"	8,4	1 x 2030 VSF	
SF 0192	DN80	10/145	2.400	1.413	1014	203	139,7	136	1/4"	10,2	1 x 3030 VSF	
SF 0432	DN100	10/145	5.040	2.966	1005	410	219,1	183	1/2"	44	3 x 2030 VSF	
SF 0576	DN100	10/145	6.720	3.955	1266	410	219,1	183	1/2"	45	3 x 3030 VSF	
SF 0768	DN150	10/145	9.600	5.650	1337	480	273	267	1/2"	70	4 x 3030 VSF	
SF 1152	DN150	10/145	13.440	7.910	1331	540	323,9	267	1"	80	6 x 3030 VSF	
SF 1536	DN200	10/145	17.200	10.124	1496	660	406,4	306	1"	200	8 x 3030 VSF	
SF 1920	DN200	10/145	21.120	12.431	1496	660	406,4	306	1"	201	10 x 3030 VSF	
											quality class - solids (ISO 8573-1)	1
											quality class - oils (ISO 8573-1)	-
											pressure drop - new element (dry) [mbar / psi]	80/1,160
											pressure drop - new element (wet) [mbar / psi]	190/2,756
											filter media	Borosilicate micro fibres
											pleated version	-
											wrapped version	✓
											sintered version	-
											min. operating temperature (°C / °F)	-20 / -4
max. operating temperature (°C / °F)	150 / 302											
CORRECTION FACTORS												
Operating pressure [bar]	2	3	4	5	6	7	8	9	10			
Operating pressure [psi]	29	44	58	72	87	100	115	130	145			
Correction factor	0.38	0.50	0.63	0.75	0.88	1	1.13	1.25	1.38			



14 bar
operating pressure

75 to 3600 Nm³/h
volume flow rate

1/4" to 3"
connections

1,5 to +150 °C
operating temperature range

stainless steel 1.4404-standard
stainless steel 1.4301-option
material

DESCRIPTION

SPF stainless steel sterile filter housings have been specifically developed for applications in process industry, where the risk for corrosion of compressed air system components is very high. To meet the required compressed air quality appropriate filter element (Sterile filter cartridge) must be installed into filter housing.

For any other technical gas please contact us or your local dealer.

SPF process filter housing can be used in variety of applications. For applications not listed please contact us or your local dealer.

For oil removal, coalescing filter element must be installed and flow direction inside-out must be provided. General arrangement is filter head on top and filter bowl on bottom.


Fluid group 1 on request.

APPLICATIONS


- Packing industry
- Biotechnology
- Breweries
- Chemical industry
- Dairies
- Fermentation processes
- Food and beverage industry
- Pharmaceutical industry
- Hospitals

SPF SERIES

STAINLESS STEEL STERILE FILTERS





TECHNICAL DATA										FILTER ELEMENTS															
Filter housing size	Pipe size D	Oper. press.	Flow rate at 7 bar(g), 20 °C		Dimensions [mm]				Mass	PI prefilter 1; 20 µm	PIW prefilter 1; 20 µm	PN prefilter 5; 25 µm	PP prefilter 3 µm	PR prefilter 1 µm	PM microfilter 0,1 µm	PS microfilter 0,01 µm	PA activated carbon	VSF sterile 0,01 µm							
	inch		bar/psi	Nm³/h	scfm	A	B	C											E	kg					
SPF 005	1/4"	14/203	75	44	225	120	76,1	1/8"	1,9	0310 PI	0310 PIW	0310 PN	0310 PP	0310 PR	0310 PM	0310 PS	0310 PA	0310 VSF							
SPF 007	3/8"	14/203	105	62	251	120	76,1	1/8"	2,0	0410 PI	0410 PIW	0410 PN	0410 PP	0410 PR	0410 PM	0410 PS	0410 PA	0410 VSF							
SPF 010	1/2"	14/203	150	88	258	121	76,1	1/8"	2,1	0420 PI	0420 PIW	0420 PN	0420 PP	0420 PR	0420 PM	0420 PS	0420 PA	0420 VSF							
SPF 018	3/4"	14/203	225	132	282	121	76,1	1/8"	2,3	0520 PI	0520 PIW	0520 PN	0520 PP	0520 PR	0520 PM	0520 PS	0520 PA	0520 VSF							
SPF 030	1"	14/203	315	185	299	136	88,9	1/8"	3,1	0525 PI	0525 PIW	0525 PN	0525 PP	0525 PR	0525 PM	0525 PS	0525 PA	0525 VSF							
SPF 047	1 1/4"	14/203	420	247	359	155	88,9	1/8"	3,4	0725 PI	0725 PIW	0725 PN	0725 PP	0725 PR	0725 PM	0725 PS	0725 PA	0725 VSF							
SPF 070	1 1/2"	14/203	600	353	395	180	114,3	1/4"	4,7	0730 PI	0730 PIW	0730 PN	0730 PP	0730 PR	0730 PM	0730-PS	0730 PA	0730 VSF							
SPF 094	2"	14/203	900	530	464	180	114,3	1/4"	5,3	1030 PI	1030 PIW	1030 PN	1030 PP	1030 PR	1030 PM	1030 PS	1030 PA	1030 VSF							
SPF 150	2"	14/203	1260	742	592	180	114,3	1/4"	6,0	1530 PI	1530 PIW	1530 PN	1530 PP	1530 PR	1530 PM	1530 PS	1530 PA	1530 VSF							
SPF 175	2 1/2"	14/203	1680	989	743	226	139,7	1/4"	11,4	2030 PI	2030 PIW	2030 PN	2030 PP	2030 PR	2030 PM	2030 PS	2030 PA	2030 VSF							
SPF 200	3"	14/203	2400	1413	995	226	139,7	1/4"	12	3030 PI	3030 PIW	3030 PN	3030 PP	3030 PR	3030 PM	3030 PS	3030 PA	3030 VSF							
SPF 240	3"	14/203	3600	2119	1029	256	168,3	1/4"	18	3050 PI	3050 PIW	3050 PN	3050 PP	3050 PR	3050 PM	3050 PS	3050 PA	3050 VSF							
										quality class - solids (ISO 8573-1)				-	-	-	6	3	2	1	1 ¹⁾	1			
										quality class - oils (ISO 8573-1)				-	-	-	-	-	2	1	1	-			
										pressure drop - new element-dry [mbar]				≤2600; ≤60		≤2600; ≤60		10	20	50	80	60	80		
										filter media				sintered stainless steel 1.4404		sintered stainless steel 1.4404		stainless steel mesh 1.4301	acrylic fibres, cellulose		borosilicate micro fibres		borosilicate micro fibres, activ. carbon		borosilicate micro fibres
										pleated version				-		-		-	✓	✓	✓	✓	-	-	
										wrapped version				-		-		✓	-	-	-	-	✓	✓	
										sintered version				✓		✓		-	-	-	-	-	-	-	
										min. operating temperature (°C / °F)				1,5 / 35		1,5 / 35		1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	1,5 / 35	-20 / -4
										max. operating temperature (°C / °F)				150 / 302		150 / 302		150 / 302	65 / 149	120 / 248	120 / 248	120 / 248	45 / 113	150 / 302	
CORRECTION FACTORS																									
Operating pressure [bar]			2	3	4	5	6	7	8	9	10	11	12	13	14										
Operating pressure [psi]			29	44	58	72	87	100	115	130	145	160	174	189	203										
Correction factor			0.38	0.50	0.63	0.75	0.88	1	1.13	1.25	1.38	1.50	1.63	1.75	1.88										

¹⁾ Valid if "S" filter cartridge is installed upstream.



AV SERIES

STAINLESS STEEL AIR VENTING FILTERS

9 to 310 Nm³/h

volume flow rate

DN32 to DN80

connections

1,5 +200 °C

operating temperature range

stainless steel **1.4404**-standard

stainless steel **1.4301**-option
material

DESCRIPTION

AV stainless steel filter housings are designed to remove impurities from air which is being feed or exhausted from tank during changing of liquid level.

To meet the required air quality appropriate filter element (typically AVF filtration grade) must be installed into filter housing.

AV filter housing is also designed for sterilisation. Before use, if needed for the application, sterilize the filters.




APPLICATIONS

- Packing industry
- Biotechnology
- Breweries
- Chemical industry
- Dairies
- Fermentation processes
- Food and beverage industry
- Pharmaceutical industry
- Water treatment systems




Sizes

Filter housing model	Connection DIN 11851	Flow rate Nm ³ /h	Dimensions [mm]		Mass kg	Filter element type
	DN		A	ø B		
AV 006	32	9	110	88,9	1,2	0310 AVF
AV 027	40	25	160	114,3	2,1	0525 AVF
AV 032	50	40	165	114,3	1,9	0530 AVF
AV 072	50	110	290	114,3	2,6	1030 AVF
AV 144	80	210	550	168,3	7,3	2030 AVF
AV 192	80	310	820	168,3	9,5	3030 AVF



Sizes - Tri-clamp connections

Filter housing model	Connection ISO 1127	Flow rate Nm ³ /h	Dimensions [mm]		Mass kg	Filter element type
			A	ø B		
AV 006	DN40 (1 1/2")	9	126	88,9	1,2	0310 AVF
AV 027	DN40 (1 1/2")	25	179	114,3	2,4	0525 AVF
AV 032	DN50 (2")	40	181	114,3	2,0	0530 AVF
AV 072	DN50 (2")	110	308	114,3	2,7	1030 AVF
AV 144	DN80 (3")	210	568	168,3	7,7	2030 AVF
AV 192	DN80 (3")	310	812	168,3	12,5	3030 AVF





1 to 3,6 bar

steam pressure

100 °C to 135 °C

steam operating temperature

up to 90 min

sterilisation time

up to 90 min

drying time

DESCRIPTION

Compressed air piping with high temperature and humidity is a suitable environment for microbial growth. In order to stop microbial activity, filters have to be sterilized. Sterilization can be done in place or in the autoclave. To perform autoclave, sterilization filters need to be removed from the system and due to septic conditions in the atmosphere it is difficult to prevent contamination during re-assembly of the installation.

Sterilization in place (SIP) ensures much better results and is more common in the process industry. Steam sterilization is an effective method to stop microbial activity. "Mobile Steam Sterilizer" is designed for easier, cheaper and faster SIP. MSS can be easily connected to any part of installation where sterilization is required. After filter sterilization is completed MSS can simply be disconnected and transported to another location.

MSS

MOBILE STEAM STERILIZER




APPLICATIONS

- Food industry
- Wine production
- Dairies
- Small breweries


INSTALLATION OF CONVENTIONAL STERILIZATION IN PLACE SYSTEM

Scheme of a conventional sterilization in place system. All the equipment surrounded by red line is required to perform sterilization of the filter housing in the upper right corner. Such a fixed steam distribution line connects steam generator to every filter.



INSTALLATION OF MSS ON STERILE FILTER

Scheme of a new Mobile Steam Sterilizer connected to a filter. After sterilization cycle is completed flexible hose for steam distribution is disconnected from the filter and complete sterilizer (equipment outlined by a red frame) can be moved to another location.



TECHNICAL DATA	MSS
Power supply	110-230V / 50/60 Hz
Power consumption	1,7 kW
Sterilization temperature	Adjustable from 100 - 135 °C
Sterilization time	Adjustable up to 90 min
Material drying time	Adjustable up to 90 min
Sound level	60 dB(A)
Dimensions W x H x -L	450 x 855 x 445 mm



Filtration grades

Below only the most common filtration grades for process applications are listed. For complete range of our products please check OMEGA AIR product catalogue or contact us for specific application.

PI grade 1; 20 μm PREFILTER particulate  Rating according to ISO 8573-1 <table><tr><td>Solid p.</td><td>Water</td><td>Oils</td></tr><tr><td>Class</td><td>Class</td><td>Class</td></tr><tr><td>-</td><td>-</td><td>-</td></tr></table> high efficient removal of solid particles from compressed air <table><tr><td>Filter media</td><td>sintered stain. steel 1.4404</td></tr><tr><td>Oper. temper.</td><td>0 - 150 °C</td></tr><tr><td>Diff. press. (dry)</td><td>60 mbar</td></tr></table>	Solid p.	Water	Oils	Class	Class	Class	-	-	-	Filter media	sintered stain. steel 1.4404	Oper. temper.	0 - 150 °C	Diff. press. (dry)	60 mbar	PIW grade 1; 20 μm PREFILTER particulate  Rating according to ISO 8573-1 <table><tr><td>Solid p.</td><td>Water</td><td>Oils</td></tr><tr><td>Class</td><td>Class</td><td>Class</td></tr><tr><td>-</td><td>-</td><td>-</td></tr></table> high efficient removal of solid particles from compressed air <table><tr><td>Filter media</td><td>sintered stain. steel 1.4404</td></tr><tr><td>Oper. temp.</td><td>0 - 150 °C</td></tr><tr><td>Diff. press. (dry)</td><td>60 mbar</td></tr></table>	Solid p.	Water	Oils	Class	Class	Class	-	-	-	Filter media	sintered stain. steel 1.4404	Oper. temp.	0 - 150 °C	Diff. press. (dry)	60 mbar	PN grade 5; 25 μm PREFILTER particulate  Rating according to ISO 8573-1 <table><tr><td>Solid p.</td><td>Water</td><td>Oils</td></tr><tr><td>Class</td><td>Class</td><td>Class</td></tr><tr><td>-</td><td>-</td><td>-</td></tr></table> high efficient removal of solid particles from compressed air <table><tr><td>Filter media</td><td>stainles steel mesh 1.4301</td></tr><tr><td>Oper. temp.</td><td>0 - 150 °C</td></tr><tr><td>Diff. press. (dry)</td><td>10 mbar</td></tr></table>	Solid p.	Water	Oils	Class	Class	Class	-	-	-	Filter media	stainles steel mesh 1.4301	Oper. temp.	0 - 150 °C	Diff. press. (dry)	10 mbar						
Solid p.	Water	Oils																																																			
Class	Class	Class																																																			
-	-	-																																																			
Filter media	sintered stain. steel 1.4404																																																				
Oper. temper.	0 - 150 °C																																																				
Diff. press. (dry)	60 mbar																																																				
Solid p.	Water	Oils																																																			
Class	Class	Class																																																			
-	-	-																																																			
Filter media	sintered stain. steel 1.4404																																																				
Oper. temp.	0 - 150 °C																																																				
Diff. press. (dry)	60 mbar																																																				
Solid p.	Water	Oils																																																			
Class	Class	Class																																																			
-	-	-																																																			
Filter media	stainles steel mesh 1.4301																																																				
Oper. temp.	0 - 150 °C																																																				
Diff. press. (dry)	10 mbar																																																				
PM grade 0,1 μm MICROFILTER particulate+coalescing  Rating according to ISO 8573-1 <table><tr><td>Solid p.</td><td>Water</td><td>Oils</td></tr><tr><td>Class</td><td>Class</td><td>Class</td></tr><tr><td>2</td><td>-</td><td>2</td></tr></table> high efficient filtration in process industry applications, removal of submicrons particles, bacteria <table><tr><td>Filter media</td><td>borosilicate micro f.</td></tr><tr><td>Oper. temperature</td><td>1,5 - 120 °C</td></tr><tr><td>Diff. pressure (dry)</td><td>50 mbar</td></tr><tr><td>Part. reten. (nom.)</td><td>99,9999% (0,1 μm)</td></tr></table>	Solid p.	Water	Oils	Class	Class	Class	2	-	2	Filter media	borosilicate micro f.	Oper. temperature	1,5 - 120 °C	Diff. pressure (dry)	50 mbar	Part. reten. (nom.)	99,9999% (0,1 μm)	PS grade 0,01 μm MICROFILTER particulate+coalescing  Rating according to ISO 8573-1 <table><tr><td>Solid p.</td><td>Water</td><td>Oils</td></tr><tr><td>Class</td><td>Class</td><td>Class</td></tr><tr><td>1</td><td>-</td><td>1</td></tr></table> high efficient filtration in process industry applications, removal of submicrons particles, bacteria <table><tr><td>Filter media</td><td>borosilicate micro f.</td></tr><tr><td>Oper. temperature</td><td>1,5 - 120 °C</td></tr><tr><td>Diff. pressure (dry)</td><td>80 mbar</td></tr><tr><td>Part. reten. (nom.)</td><td>99,9999% (0,01 μm)</td></tr></table>	Solid p.	Water	Oils	Class	Class	Class	1	-	1	Filter media	borosilicate micro f.	Oper. temperature	1,5 - 120 °C	Diff. pressure (dry)	80 mbar	Part. reten. (nom.)	99,9999% (0,01 μm)	PA grade ADSORPTION activated carbon  Rating according to ISO 8573-1 <table><tr><td>Solid p.</td><td>Water</td><td>Oils</td></tr><tr><td>Class</td><td>Class</td><td>Class</td></tr><tr><td>-</td><td>-</td><td>-</td></tr></table> high efficient removal of oil, hydrocarbons, vapours and odours from compressed air <table><tr><td>Filter media</td><td>activated carbon, borosilicate micro f.</td></tr><tr><td>Oper. temperature</td><td>1,5 - 45 °C</td></tr><tr><td>Diff. pressure (dry)</td><td>60 mbar</td></tr><tr><td>Residual oil cont.</td><td><0,0005 mg/m³</td></tr></table>	Solid p.	Water	Oils	Class	Class	Class	-	-	-	Filter media	activated carbon, borosilicate micro f.	Oper. temperature	1,5 - 45 °C	Diff. pressure (dry)	60 mbar	Residual oil cont.	<0,0005 mg/m³
Solid p.	Water	Oils																																																			
Class	Class	Class																																																			
2	-	2																																																			
Filter media	borosilicate micro f.																																																				
Oper. temperature	1,5 - 120 °C																																																				
Diff. pressure (dry)	50 mbar																																																				
Part. reten. (nom.)	99,9999% (0,1 μm)																																																				
Solid p.	Water	Oils																																																			
Class	Class	Class																																																			
1	-	1																																																			
Filter media	borosilicate micro f.																																																				
Oper. temperature	1,5 - 120 °C																																																				
Diff. pressure (dry)	80 mbar																																																				
Part. reten. (nom.)	99,9999% (0,01 μm)																																																				
Solid p.	Water	Oils																																																			
Class	Class	Class																																																			
-	-	-																																																			
Filter media	activated carbon, borosilicate micro f.																																																				
Oper. temperature	1,5 - 45 °C																																																				
Diff. pressure (dry)	60 mbar																																																				
Residual oil cont.	<0,0005 mg/m³																																																				

PP grade**3 µm****PREFILTER**
particulate**Rating according to ISO 8573-1**

Solid p.	Water	Oils
Class	Class	Class
6	-	-

high efficient removal of solid particles from compressed air

Filter media	acrole fibres, cellulose
Oper. temperature	1,5 - 65 °C
Diff. press. (dry)	10 mbar
Part. reten. (nom.)	99,99% (3 µm)
Part. retent. (ISO)	95 %

PR grade**1 µm****PREFILTER**
particulate**Rating according to ISO 8573-1**

Solid p.	Water	Oils
Class	Class	Class
3	-	-

high efficient removal of solid particles from compressed air

Filter media	borosilicate micro f.
Oper. temperature	1,5 - 120 °C
Diff. pressure (dry)	20 mbar
Part. reten. (nom.)	99,9999% (1 µm)
Part. retent. (ISO)	99.8 %

VSF grade**STERILE**
particulate + bacteria removal +
sterile**Rating according to ISO 8573-1**

Solid p.	Water	Oils
Class	Class	Class
1	-	-

high efficient sterile filtration of compressed air, process air and technical gasses

Filter media	borosilicate micro f.
Oper. temperature	-20 to +150 °C
Diff. pressure (dry)	80 mbar
Part. reten. (nom.)	99,9999% (0,1 µm)

AVF grade**STERILE**
particulate + bacteria removal +
sterile**Rating according to ISO 8573-1**

Solid p.	Water	Oils
Class	Class	Class
1	-	-

high efficient sterile filtration of compressed air, process air and technical gasses

Filter media	borosilicate micro f.
Oper. temperature	-20 to +150 °C
Diff. pressure (dry)	80 mbar
Part. reten. (nom.)	99,9999% (0,1 µm)

OMEGA AIR

Air and Gas



OMEGA AIR d.o.o. Ljubljana

T +386 (0)1 200 68 00
F +386 (0)1 200 68 50

info@omega-air.si

Cesta Dolomitskega odreda 10
SI-1000 Ljubljana, Slovenia
www.omega-air.si

950119 - 10/2019

